# monolith

SISTEMA DE IMPLANTES  
IMPLANT SYSTEM

[www.gmidental.com](http://www.gmidental.com)

Implantes dentales de conexión dual  
Dual connection dental implants

Made in Spain

 gmi

## Índice - Summary

### CARACTERÍSTICAS GENERALES

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### KIT QUIRÚRGICO GMI MONOLITH

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### KEY FEATURES

- Surface treatment / studies 4
- Recommendations of use 6
- Labelling description 7

### IMPLANT SYSTEM

- **gmi monolith** dental implants 8
- Insertion protocol for **gmi monolith** implants 9
- **gmi monolith** implants range 10
- Drilling sequences for **gmi monolith** implants 11

### ATTACHMENTS

- **gmi monolith** healing abutments 12
- **gmi monolith** carrier abutments 13
- **gmi monolith** ball abutments 14
- **gmi monolith** angled abutments 15
- **gmi monolith** straight abutments 15
- **gmi monolith** replicas 16
- **gmi monolith** impression coping 16
- **gmi monolith** clinic and laboratory screws 17
- Recommended tightening torques 18

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


To promote the adhesion and growth of bone cells, the outer surface of the implant has been treated with the exclusive ADS® treatment (Advanced Double-Grip Surface) which combines a white corundum micro-bubble treatment and acid etching, for a non-smooth roughness that maximizes the contact surface between implant and bone and therefore improves stability and adhesion between them.

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
The ADS® treatment is the result of a series of collaborations with institutes and universities, as well as our own research and tests, that allow us to monitor and validate the optimal biological response to the implant surface: chemical (X-ray scattering studies to find out the composition of the implant and its surface) and topography studies (roughness studies) have been carried out, as well as biological (in vivo response studies in animals) and clinical studies.

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X-ray scattering studies show that the surface of the implants is only made of titanium and oxygen, which are the basic components of the oxide layer that is formed in the implant passivation process and that gives it its excellent corrosion resistance.

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
In order to control and identify the topography of the implant surface, roughness studies were performed with profilometer, as well as with techniques such as SEM (Scanning Electron Microscope) and CLSM (Confocal Laser Scanning Microscope), allowing to obtain and control the mean roughness (Ra) values within the parameters recommended in international publications.




Perfil de rugosidad obtenido con HOMMEL-WERKE T1000 • Roughness profile obtained with HOMMEL-WERKE T1000.

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To determine the biological response of **gmi** implants, in vivo studies were conducted by inserting the implants on animals, leaving them without load during the healing process and performing a histological study to prove the excellent response of the bone cells and the complete osseointegration of the implant.



SEM-BSE: detalles ultraestructurales de las células en biofilm sobre el implante. • SEM-BSE: Ultrastructural details of cells in biofilm on the implant.

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The **gmi monolith** implant system was designed and tested to be used through one-stage surgery in upper or lower jaw in cases of total or partial edentulism or in cases of loss of isolated teeth (single restoration), with a 12-degree maximum inclination to the perpendicular plane to the occlusal plane, in accordance with the following diagram:

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

The **gmi monolith** implant system allows performing the following restorations:

- Unitary cemented restorations of mandibular or maxillary lateral incisor. For these cases, h=2 mm implant must be used.
- Multiple cemented restorations of small bridges with splinted implants in incisor. For these cases, h=2 mm implant must be used.
- Stabilization of implant-retained and tissue-supported overdentures. For these cases, h=2 mm and h=4 mm implants can be used.

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
The **gmi monolith** implant range is designed to be inserted applying a 35-40 N-cm torque, although it can bear a maximum tightening torque of 60 N-cm. **gmi** recommends manual insertion of the implant with a dynamometric ratchet wrench to ensure that these values are not exceeded.

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The **gmi monolith** range of implants and prosthetic attachments are designed and manufactured with tolerances intended to maintain adjustments that minimize the clearance between components and therefore guarantee the stability of the connection. **gmi** recommends the use of original attachments to ensure a perfect fit between the components of the restoration.

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

To prevent prior deformation of the clinic screw, **gmi** recommends restricting its use only for the final tightening of the prosthesis. To make adjustments during the prosthetic procedures, use of the laboratory screw is recommended.

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
In the chart are detailed all the symbols appearing on the **gmi monolith** implant system's labelling and packaging and their corresponding description.

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\* **Minimum requirements for downloading and opening the electronic instructions for use:**  
 - Stable internet connectivity.  
 - A device with a modern operating system that supports PDF viewing.  
 - Adequate storage space on your device.  
 - A web browser such as Google Chrome, Mozilla Firefox, Safari, Microsoft Edge or Opera installed on your device.  
 - A PDF reader application installed on your device.

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
gmi **monolith** dual connection implants are machined in c.p. grade IV titanium of the highest quality, thus ensuring a high chemical stability and excellent biocompatibility. The titanium oxide layer, that is formed on the implant surface after the passivation process, and the careful design of the outer surface of the thread, guarantee good osseointegration and high primary stability.

The external dual inlet thread of the implant helps to reduce surgery time. Moreover, self-threading millings and slight tapering facilitate implant insertion, reducing tension at the bone-implant interface and thus preventing threading stress-related problems.

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⑧ Rosca cuerpo - 2 entradas  
Paso: 1,60 mm

① External square  
B/F 2.3 mm

② Dual connection  
External cone of 5°  
Internal cone of 5°

③ Gingival heights  
2 / 4 mm

④ Length  
12 / 14 / 16\* mm

⑤ Self threading millings  
2 millings at 180°

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**gmi monolith** dual connection implants are packed in a screen-printed cardboard box, with 2 traceability labels and a double packaging containing the implant system, with the secondary packaging made of technical plastic and the primary packaging made of Pyrex glass to ensure the inert performance of the container. Once packaged in a clean room, they are sterilized with gamma radiation according to regulations.

The packaging system of the **gmi monolith** implants consists of a titanium support which houses the implant-implant carrier set, prepared to be inserted directly into the mouth, using a single hex-2.50 mm implant carrier wrench, designed for this purpose. This new system offers the odontologist the following advantages: more control and a better view of the implant insertion process, it is easier to use in reduced interdental spaces and a reduction in surgery time.\*

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\* See the instructions for use in electronic format.


1. Check sticker indicating implant length and open the tab of the carton.
2. Remove the secondary packaging and product identification labels.
3. Check the integrity of the safety seal and unscrew the plug of the secondary packaging.
4. Remove the primary packaging from the secondary packaging avoiding hitting it against a hard surface.
5. Extract primary plug making a lateral movement.
6. Insert the implant carrier wrench until a slight retention is noticed, addressing the faces of the key with the faces of the implant carrier hex broaching.
7. Check that the key is fully inserted and turn slightly while gently pulling on the assembly.
8. Place the implant in the bone bed by applying a torque of 35 N-cm until the treated area is at crestal or subcrestal level. Remove the key.
9. Untighten screw using HEX-1.20 mm wrench and remove it from the implant carrier.
10. Introduce the extraction wrench and remove the implant carrier.

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

The following is a description of the drilling sequences and the recommended drilling conditions for their use:

- Lance-shaped: 1,200 - 1,500 rpm.
- Step drills:
  - Ø2.00 / Ø2.50 mm → 900 - 1,200 rpm.
  - Ø2.50 / Ø2.80 mm → 500 - 700 rpm.


• For drilling, external cooling with saline solution is required.

monolith Ø3,00 - L12		
	Descr.	Ref.
1	Fresa lanza	Lance-shaped drill
2	Medidor de encía	Soft tissue gauge
3	Fresa escalonada Ø2,00 / Ø2,50	Step drill Ø2.00 / Ø2.50
4	Fresa escalonada Ø2,50 / Ø2,80	Step drill Ø2.50 / Ø2.80

monolith Ø3,00 - L14		
	Descr.	Ref.
1	Fresa lanza	Lance-shaped drill
2	Medidor de encía	Soft tissue gauge
3	Fresa escalonada Ø2,00 / Ø2,50	Step drill Ø2.00 / Ø2.50
4	Fresa escalonada Ø2,50 / Ø2,80	Step drill Ø2.50 / Ø2.80

monolith Ø3,00 - L16		
	Descr.	Ref.
1	Fresa lanza	Lance-shaped drill
2	Medidor de encía	Soft tissue gauge
3	Fresa escalonada Ø2,00 / Ø2,50	Step drill Ø2.00 / Ø2.50
4	Fresa escalonada Ø2,50 / Ø2,80	Step drill Ø2.50 / Ø2.80

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The healing abutment is threaded onto the implant and it is responsible for preparing a mucous way for connecting the implant to the secondary structure or prosthesis, once osseointegration is complete. The healing abutment should stay in place for a minimum of two to three weeks until the wound has finished healing.

### Key features

- Machined body in high strength PEEK.
- Machined grade V titanium screw according to ISO-5832.
- Tightening with b/f 1.20 mm hex wrench.
- Available in emergency height of 5 mm.
- Manual tightening torque: Max. 15 N-cm.
- Self-unlocking system (auto-extractor).

Ø Cabeza (mm) Ø Head (mm)	Altura (mm) Height (mm)	Ref.
Ø 4,00	5,00	KDPoF5001
Ø 5,00	5,00	KDPoF5002

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The carrier abutment is defined as a machined attachment that once fixed with the clinical screw to the implant, serve as a support to make cemented restorations. As a complementary function, carrier abutments allow to make the closed tray impression coping using the suitable plastic abutment.

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### Key features

- Grade V titanium machining as per ISO-5832 standard.
- Compatible with b/f 2.50 mm hex wrench to use as an implant carrier.
- Clinic screw included.
- Screw tightening torque: 25 N-cm.

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Ball abutment is defined as a machined attachment being fixed to the top of the implant and serving as the retention element for a removable prosthesis, with the previous addition of the retention caps that will be fixed to the top spherical part of the ball abutment.

### Key features

- Grade V titanium machining as per ISO-5832 standard.
- Yellow extra-soft retention cap and metal container included.
- Protection disk included.
- Tightening with b/f 1.20 mm hex wrench.
- Instruments for caps: Insertion / removal tool for caps (RH485IC).
- Tightening torque: 30 N-cm.

Ø Bola (mm) Ø Ball (mm)	Ref.
Ø 2,50	KDFoF5001

Retenciones Retention caps			
Tipo retención Retainer type	Color	Cantidad (ud.) Quantity (Un)	Ref.
Extra-suave • Extra-soft (500 gr)	Amarillo • Yellow	6	RH60CRNAY
Suave • Soft (900 gr)	Rosa • Pink	6	RHCN001
Estándar • Standard (1300 gr)	Transparente • Transparent	6	RHCN002
Elástica • Elastic (350 gr)	Verde • Green	6	RHCN031

Contenedor cofias Metal container	
Cantidad (ud.) Quantity (Un)	Ref.
2	RH041CANB

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When it is necessary to correct an extreme angular position of the implant inside the bone structure in relation to adjacent natural or prosthetic teeth, angled abutments should be used. Once mounted and oriented on the implant, these abutments will allow the parallelization of the new crown with other structures.

### Key features

- Grade V titanium machining as per ISO-5832 standard.
- Allow non-parallelism correction up to 12 degrees.
- Clinic screw included.
- Abutment shoulder adapted to the shape of the gum.
- Screw tightening torque: 25 N-cm.

Ang. (°)	Altura (mm) Height (mm)	Ref.
12°	9,20	KDMoF5002

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Once the mucosal tunnel has been formed by the healing abutment and after removing it from this location, it is appropriate to use the abutment, which is defined as the machined tapered attachment that, after being fixed with the clinic screw to the top of the implant, acts as a support for the clinic crown and is attached to it using cement.

### Key features

- Grade V titanium machining as per ISO-5832 standard.
- Anti-rotational planes to increase retention.
- Available only in non-rotary version.
- Clinic screw included.
- Screw tightening torque: 25 N-cm.

Ø Hombro Ø Shoulder	Altura (mm) Height (mm)	Ref.
Ø 3,80	9,00	KDMoF5001

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Once the transfer of the position has been made using the impression coping from the biological medium to the laboratory working model, it is necessary to use replicas, which are defined as the machined attachments designed on one hand to reproduce the implant position in the mouth of the patient to a working model, and on the other to serve as a connection model for the construction in the laboratory of the structure to carry out the restoration.

### Key features

- Grade V titanium machining as per ISO-5832 standard.
- Anti-rotary millings and vertical retention neck for a good fixation.

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Impression coping are defined as the machined attachments that are used to make the transfer of the position of the implant in the biological medium to a laboratory model, on which the prosthetic will work to create the dental crowns or prostheses. This process is achieved using impression materials placed in an appropriate manner, that will harden in the oral cavity and produce the negative of the patient's mouth. Once extracted, replicas are placed and the casting plaster is poured to obtain the positive mould where the replica will be placed in the original position of the implant in the mouth.

### Key features

- CT impression coping machined in technical plastic.

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Clinic screws are defined as machined devices used for the final anchoring of the abutments to the implant, while laboratory screws are the machined devices used for the temporary anchoring of abutments to the replica of the implant used by the dental laboratory where the construction of dental crowns or prostheses takes place.

### Key features

- Machined grade V titanium screws according to ISO-5832.
- Tightening of screws with b/f 1.20 mm hex wrench.
- Tightening torque: 25 N-cm.

Tornillo clínica  
Clinic screw

Hex. llave (mm)  
Hex. wrench (mm)

Ref.

⊙1,20 mm

KDT0C5001

Tornillo laboratorio  
Laboratory screw

Hex. llave (mm)  
Hex. wrench (mm)


Ref.

⊙1,20 mm

KDT0C5002

## Pares recomendados de apriete - Recommended tightening torques

Par máximo inserción implante • Implant insertion maximum torque

Apriete tornillo clínica • Clinic screw tightening torque

Apriete pilar bola • Ball abutment tightening torque

Apriete pilares de cicatrización • Healing abutments tightening torque

## Kit quirúrgico gmi monolith - gmi monolith surgical kit

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**gmi monolith surgical kit** consists of an autoclavable case made of high strength resin, containing all the necessary material for the correct positioning of **gmi monolith** dual connection dental implants system. The components are arranged in such a way to be easy to find in each moment and in any stage of the intervention.

**gmi surgical kits** are distributed in cartons, made of high strength technical polymer, and should be steam autoclaved before use at a maximum temperature of 134°C.

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